




Revised 2019

Prep device for installation:

Before starting install, determine guide locations along the edge of the opening. Guide locations should be evenly spaced (usually every other hole).

Mark and punch holes.



Install Hints:

Use QuickFit™ Tool 1/8" punch to make required holes.

If using existing holes for rivet guides, remove metal or plastic grommets from eyelets to ensure a secure rivet setting.

Rivets can be hand snapped to position for verifying fit, but system must be set with QuickFit™ Tool.

Press rivets at 90° with QuickFit™ Tool; pressing an angle will result in poor rivet setting.

Secure dial on thicker/sturdier part of tongue.

Recommended Tools:

- QuickFit™ Setter Tool
- QuickFit™ Rivet Pack
- Marker Pen




Advisory

Must be installed by a trained professional.
Lace will wear. At minimum, inspect monthly.
Overtightening could reduce circulation.
Do not use on patients with nerve issues.
Keep lace away from sharp edges.


See advisories in other languages: clickmedical.co/advisories

	WARRANTY US customers contact Click Medical directly: Help@ClickMedical.co Help Line: +1.970.670.7012 International customers, please contact your local distributor	<table border="1"> <tr> <td data-bbox="1157 936 1252 989"> EC REP </td><td data-bbox="1252 926 1386 999"> MDSS GmbH Schiffgraben 41 30175 Hannover, Germany </td><td data-bbox="1386 936 1511 989"> </td></tr> <tr> <td data-bbox="1157 999 1252 1100"> </td><td colspan="2" data-bbox="1252 999 1511 1100"> Click Medical, LLC 1205 Hilltop Parkway, W101 Steamboat Springs, CO, 80487, USA +1-970-670-7012 </td></tr> </table>	EC REP	MDSS GmbH Schiffgraben 41 30175 Hannover, Germany			Click Medical, LLC 1205 Hilltop Parkway, W101 Steamboat Springs, CO, 80487, USA +1-970-670-7012	
EC REP	MDSS GmbH Schiffgraben 41 30175 Hannover, Germany							
	Click Medical, LLC 1205 Hilltop Parkway, W101 Steamboat Springs, CO, 80487, USA +1-970-670-7012							

- 1 Peel sticker from its backing paper. Use sticker to position the system on the device. Allow all the guides to freely move.




- 2
 - a) Attach lace guides by inserting male rivet from inside the device and through lace guide.
 - b) Snap female rivet cap on top of lace guide, holding the guide to the shoe.




Hint: Start with guides farthest away from the dial. Repeat for all guides.


- 3 Remove sticker to determine placement of dial. Allow 1cm of clearance between top of tongue and top of dial. Mark location of holes.



- 4 Punch hole using the 1/8" punch option of QuickFit™ Tool.




- 5 Snap rivets as described in step 2. If satisfied with placement, secure using QuickFit™ Tool.



- 6 Repeat step 5 for all guides.



- 7 Check alignment of guides.



- 8 Cycle system to test for smooth function:
Push to engage.
Twist to tighten.
Pull to release.

